



SORT™

A standardized approach to IV line management

Four steps. Every line. Every time.

Organize tubing from pump to patient using one consistent workflow.

START CLEAN Perform hand hygiene • Ensure adequate lighting • Visualize the entire IV pathway from pump to patient

S

SEPARATE



- Insert one IV line per organizer channel
- Keep lines parallel with no crossing or tangling
- Maintain slack to prevent tubing tension

HELPS REDUCE LINE CROSSOVER

O

ORGANIZE



- Use a minimum of two LINE-SORT pads per set of IV lines, approximately 12–15 inches apart along the tubing pathway, depending on the clinical setup
- Label the IV bundle per facility policy
- Reassess placement after repositioning, transfer, or transport

STANDARDIZES TUBING ORGANIZATION

R

RECONCILE



- Confirm every line matches its active medication or fluid order, per facility policy
 - ✓ Medication/solution
 - ✓ Pump/channel
 - ✓ Route
 - ✓ Connection site
 - ✓ Patient
- Mismatch? Stop, correct, recheck

SUPPORTS ORGANIZED LINE IDENTIFICATION

T

TRACE



- Follow your facility's tracing protocol
- Trace and verify IV lines before administration or adjustment, in accordance with facility policy
- Hand-over-hand, **Patient → Source**
- Confirm the correct line before use

SUPPORTS EASIER LINE TRACING

RUN SORT during line setup, line additions, new infusions, transfers, shift changes, handoffs, and after tubing reconfiguration.

Always follow your facility's policies and procedures. LINE-SORT is an organizational aid and does not replace clinical judgment, institutional protocols, or standard IV line tracing procedures. See the Directions for Use for complete instructions, warnings, and precautions.

Clarity in seconds. When minutes matter most.